

THE AIR JET BODY MASSAGE DEVICE DEDICATED TO BEAUTY AND HEALTH



Aeroline® is a device for professional use in esthetical and medical massage therapies.

The basic working principle of the Aeroline® device is based on a massage action of a strong air flow, which is directed by a masseur i.e. by a therapist to the part of the body which is being treated by a massage.

By a professional directing of the air flow the massage effects are achieved which cannot be achieved by any of the so far known massage technologies.

The basic advantages of this new body massage technology by air flow – air massage are as follows;

- A totally natural method,
- A very comfortable massage,
- Achievement of extraordinary results,
- Use in various sorts of therapies,
- Speed – the massage treatment of the whole body lasts for only 25 minutes,
- Multiple reduction in physical fatigue of the masseur–therapist,
- During the massage treatments there is no direct physical contact between the masseur – therapist and the treated person.

The most frequent use of the Aeroline® air massage:

- Reduction of body weight,
- Anti-cellulite treatments,
- Body sculpturing,
- Firming of the muscular mass,
- Lymphatic drainage,
- Relaxing treatments,
- Poor periphery circulation,
- Rheumatism,
- Treatments following a post-trauma immobilization,
- Sports injuries.

The air massages by the Aeroline® device are most frequently performed in a set of 10 treatments over a two – weeks period.

The wanted results are usually achieved already after one or two massage sets.

The number of massage treatments by means of the air flow is not limited, so that the massage can be applied until achievement of the wanted results.

In order to maintain the achieved condition one massage per week is usually sufficient.

The therapies by means of the Aeroline® device, i.e. the massage by the air flow, can be performed only by an educated cosmetician, masseur or physiotherapist.

Adequate working room for locating the Aeroline® air massage is a smaller room or cabin with minimal dimensions of 2,2 x 2,8 meters.

The working room must have an inflow of fresh air – a smaller window, or a high quality ventilation system.

The Aeroline® device has been constructed in such a way not to occupy the working space but it is situated underneath the massage bed on which the treatments are being performed. Most usually, a classical internist bed is used, which has a free room underneath the working space of minimum 62 cm.



Guarantee: 2 years

Maintenance: Regular service every 12 months / 1000 work hours.

Technical data	
Flow pressure	0,86 N/cm ²
Area of air flow	11 cm ²
Temperature of air	38 – 45 °C
Voltage	200-240 V 50/60Hz
Power	< 2,5 kW
Weight	36 kg
Dimensions	60x70x61 cm

PCT WO/2007/003974

The producer reserves the rights of change in technical characteristic



PRINTED IN CROATIA

CROAL SMOK
d.o.o. za proizvodnju trgovinu i usluge

Novomarska 32
10000 Zagreb CROATIA
E-mail: info@aeroline.com.hr
Website: www.aeroline.com.hr